

Apex BI Product List

Mesa Labs Vaporized Hydrogen Peroxide Biological Indicators

APEX Biological Indicator products are designed specifically for the pharmaceutical, food and medical device industries utilizing H₂O₂ vapor sterilization. The stainless steel carrier is designed for use with today's isolator and filling line applications.

Item	Reorder Number	Population	ATCC Number	Package Configuration
Stainless Steel Discs in Tyvek/Tyvek <i>G. stearothermophilus</i>	HMV-091E3	10x3	12980	100/pkg
	HMV-091	10x6	12980	100/pkg
	KCD-404	10x4	12980	100/pkg
Bare Metal BI Stainless Steel Ribbon <i>G. stearothermophilus</i>	SBC-327	10x6	12980	100/pkg
Tri-Scale BI; Tyvek pkg BI w/ 10x6, 10x5, 10x4 Stainless Steel Carriers <i>G. stearothermophilus</i>	LOG-456	10x4, 10x5, 10x6	12980	50/pkg
Stainless Steel Discs in Tyvek/Tyvek <i>B. atrophaeus</i>	GRS-090E3	10x3	9372	100/pkg
	GRS-090	10x6	9372	100/pkg
Stainless Steel Discs in Tyvek/Tyvek <i>G. stearothermophilus</i>	NAS-152E3	10x3	7953	100/pkg
	NAS-152	10x6	7953	100/pkg
Sterile Tyvek Transfer Pockets for SBC-327	BMWK-750			100/pkg
Spore Suspension for gaseous hydrogen peroxide <i>G. stearothermophilus</i>	LAU-196	10x7/0.1 ml	12980	10 ml vial
	LPT-606	10x7/0.1 ml	7953	10 ml vial
Releasat Purple, for use with HMV-091, KCD-404, LOG-456, NAS-152, SBC-327 3.6 Fill Volume in 16mm x 100mm Tubes	PM/100	--	--	100 units

For more information, contact Mesa Labs or your local distributor.

