

Product Specification

Releasat® Purple Media

Product Description

Releasat® Purple media is specially formulated soybean casein digest culture medium containing the pH indicator Bromocresol Purple. The Releasat culture medium is supplied in individual tubes for easy aseptic transfer and culturing. Test results are obtained during incubation by observing the medium which changes color from purple to or toward yellow and/or turbid if any spores have survived the process.

Technical Parameters

Incubation	Time: 7 days
Storage	Room temperature. Protect from light.
Shelf-Life	12 months from date of manufacture
Nominal Dimensions	Media Tube: 16 x 100mm Volume: 3.6mL

Indications for Use

Releasat Purple media is specially formulated for rapid outgrowth of *Geobacillus stearothermophilus* spores that may have survived steam and hydrogen peroxide sterilization/decontamination processes.

Packaging

Releasat Purple media is packaged as 100 tubes of Releasat culture media in a chipboard box.

Labeling

- Each media tube is labeled with lot number and expiration date.
- Each box has a self-adhesive label on the outside to further identify the product.

Reorder Information

- PM/100 (100 sets/box)

